**INSTRUCTIONS FOR USE****Description:**

The Emboline Angiographic Catheter is a single-lumen catheter that consists of a winged hub, strain relief, and radiopaque catheter shaft. Angiographic catheters are single-use, sterile, non-pyrogenic, intravascular catheters intended for short-term use.

**User(s):** For use by physicians trained in diagnostic and interventional radiology, cardiology, nephrology, and vascular surgery procedures.

**Patient Population:**

Emboline Angiographic Catheters are intended for use in adult patients undergoing intravascular diagnostic procedures by trained physicians.

**How Supplied:**

Emboline catheters are packaged in a straight position with a flush port to facilitate compliance with the manufacturer recommended guidelines that the catheter be flushed with saline or heparinized saline prior to use (See directions for use - Note). Individually packaged 5-10 per box, refer to catalog for ordering information.

**Intended Purpose:**

The Emboline Angiographic Catheter is intended to be used to deliver radiopaque media to selected sites in the vascular system.

**Indications for Use:**

The Emboline Angiographic Catheter is indicated for use in patients undergoing cardiovascular angiography procedures. Vascular access insertion sites are pre-determined by the physician performing the procedure based on the patient vascular disease process and anatomical vessel size.

**Contraindications:**

The Emboline Angiographic Catheters have no known contraindications.

**CLINICAL BENEFITS:**

Emboline Angiographic Catheters provide an indirect clinical benefit to patients by facilitating delivery of contrast media during routine diagnostic procedures. The device has a clinical lifetime of use for the duration of the associated procedure, during which the clinical benefit is supplied for as long as the device is in use.

**PERFORMANCE CHARACTERISTICS:**

Emboline Angiographic Catheters are designed with performance characteristics for use in a patient's vascular system. See individual labels for French size, Length, Flow rate, Maximum guide wire size, and Pressure ratings.

**Warnings:**

- Angiographic catheters are single use.
- Recommended guide wire size: See the individual catheter label.

- Angiographic catheters may be used with syringes or power injectors up to the labeled flow rate and pressure rating.
- Guide wire use has been associated with greater incidence of thrombus formation.
- Never withdraw the catheter against resistance until the cause of the resistance is determined by fluoroscopy. Movement of the catheter against resistance may result in catheter damage leading to potential catheter tip separation.
- **Pigtail Catheters.** Ensure the pigtail straightener is removed from the catheter shaft by gently pulling on the tear-away tab handle prior to advancing the pigtail into the introducer sheath.
- Anticoagulation therapy, per facility protocol, should be considered to reduce potential for thrombus formation on the device.
- In the event of a malfunction of the device and/or changes in the performance of the device, exercise caution as this may indicate a change that may affect the safety of the device.
- After use, dispose of device in a manner consistent with standard protocols for biohazard waste disposal.

**Precautions:**

- Inspect all catheters prior to use. Do not use if package is opened or damaged.
- Do not use after expiration date.
- Do not autoclave.
- Confirm compatibility of the catheter and other associated devices to be used prior to use.
- Angiography should be undertaken only by an experienced angiographer. Physicians should be familiar with the use of Angiographic products and the literature concerning the complications of angiography.
- For high pressure contrast delivery do not exceed maximum PSI ratings.
- Fluids intended for injection per the indications for use are radiographic contrast, heparinized saline, and saline. Do not inject other fluids without confirmation of compatibility from the fluid manufacturer.
- If resistance is felt when removing the guidewire from the catheter, remove the guide wire and catheter as a unit to prevent potential damage to the vessel wall.
- Exposure to VHP vaporized Hydrogen Peroxide may cause catheter failure.
- In the EU, any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the applicable Member State.

**Potential Complications:**

Emboline Angiographic Catheters should be used by physicians familiar with percutaneous catheter introduction. Complications which may be associated with the use of angiographic catheters include, but are not limited to, the following:

- Thrombus formation/emboli
- Plaque dislodgement
- Vascular occlusion
- Stroke
- Death
- Air embolism
- Arterial wall damage
- Infection
- Vascular spasm
- Foreign Body / Tip Fracture

Some of the stated potential adverse events may require additional surgical intervention.

Directions for Use:

1. A coated guide wire is recommended for use with Emboline Angiographic Catheters.
2. Utilize aseptic technique, carefully remove the catheter from the package by grasping the hub and slowly withdrawing it from the package.
3. Flush all devices entering a blood vessel with sterile heparinized saline or similar isotonic solution before each use.
4. Introduce the catheter into the blood vessel using vascular entry technique of choice.
5. Do not use excessive force to advance or manipulate a catheter through the vessel. Advancement with excessive force may cause vascular damage with the associated device (guide wire).
6. If resistance is felt during advancement, stop immediately to identify the source of resistance, and confirm the tip position under fluoroscopy.
7. Refer to package label for flow rates and maximum pressure ratings.
8. Remove the catheter and guide wire as a single unit to avoid vessel wall damage.
9. **Pigtail Catheters:** Pigtail catheters are fitted with a curve straightener. It is split for easy removal after the catheter tip is on the guide wire.
  - a. Slide straightener forward until pigtail is straight.
  - b. Place catheter tip on guide wire and remove straightener by pulling on the straightener tab.

Sterile and non-pyrogenic in unopened and undamaged inner package. Store in a cool dry place.

REUSE PRECAUTION STATEMENT

For single patient use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness, or death. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness, or death of the patient.

For a copy of this device's current European Summary of Safety and Clinical Performance (SSCP), please go to the European database on medical devices (Eudamed), where it is linked to the basic UDI-DI. <https://ec.europa.eu/tools/eudamed>.

Basic UDI-DI: 08600057436EPT-01QK

| SYMBOL                                                                             | DESIGNATION                                                                                                                                                                                                                  |
|------------------------------------------------------------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
|    | Catalog Number                                                                                                                                                                                                               |
|   | Lot Number                                                                                                                                                                                                                   |
|   | Caution                                                                                                                                                                                                                      |
|   | Consult Instructions for Use. For electronic copy go to <a href="http://www.emboline.com/contact-us">www.emboline.com/contact-us</a> . For printed copy, contact Emboline via the contact information provided on the labels |
|   | Sterilized using Ethylene Oxide                                                                                                                                                                                              |
|   | Caution: Federal Law (USA) restricts this device to sale by or on the order of a physician.                                                                                                                                  |
|   | Do not use if package is damaged.                                                                                                                                                                                            |
|   | Single Use                                                                                                                                                                                                                   |
|   | Use By                                                                                                                                                                                                                       |
|   | Do Not Resterilize                                                                                                                                                                                                           |
|   | Non-pyrogenic                                                                                                                                                                                                                |
|   | Single sterile barrier system                                                                                                                                                                                                |
|   | Manufacturer                                                                                                                                                                                                                 |
|   | Medical Device                                                                                                                                                                                                               |
|   | Date of Manufacture                                                                                                                                                                                                          |
|   | Unique Device Identifier                                                                                                                                                                                                     |
|  | Authorized Representative in European Community                                                                                                                                                                              |

|                                                                                             |                                                       |
|---------------------------------------------------------------------------------------------|-------------------------------------------------------|
| <div>F</div>                                                                                | French Size                                           |
| <div>PSI</div>                                                                              | Recom.pressure Limit                                  |
| <div>GW</div>                                                                               | Recom.guidewire                                       |
| <div>SH</div>                                                                               | Side Holes                                            |
| <div>FLO</div>                                                                              | Rates Achieved at 15.0cP individual Results May Vary. |
| <div>QTY</div>                                                                              | Quantity                                              |
| <div>CM</div>                                                                               | Length                                                |
| <div></div> | Keep Dry                                              |
| <div></div> | Keep Away from Sunlight                               |

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Emboline, Inc.
AR Experts B.V.